

COMBINATION FREEZER & REFRIGERATOR

MODEL CFR-325



Model CFR-325

Features and Benefits

- CFC free refrigeration – environmental friendly.
- Built-in overload protection for compressor – ensures lasting operation.
- Lock provided as standard – secure against unauthorized access.
- Door with fan switch stops circulation fan when door is open-reduces temperature fluctuations.
- Large stainless steel cabinet – provides extra storage space and durable usage.
- Temperature control – factory set at -25°C (adjustable from -15°C to -25°C for Freezer).
- Temperature control – factory set at 4°C (adjustable from 2°C to 8°C for Refrigerator).
- Digital thermometer – temperature is display numerically.
- Automatic condensate removal – eliminates the need for drain lines.
- Forced air circulation – keeps temperature consistent at all levels.
- Double doors – separate freezer and refrigerator compartments.
- Runs on demand refrigeration design – reduces power consumption.
- Full – length adjustable shelves – for easy loading and unloading of stored products.
- Caster wheel with locking brakes – for easy transportation and safety at site when placed in position.

Specifications

MODEL	CFR-325	
Type	Combination - freezer and refrigerator	
Cabinet construction	Inner wall and outer casing of stainless steel	
Cabinet insulation	CFC-free polyurethane, foam-in-place insulation	
Temp. working range	-15°C to -25°C (Freezer), 2°C to 8°C (Refrigerator)	
Freezer	Controller	Digital microprocessor controller with defog, lock and off-set functions, Hi-Low alarm
	Temp. setting & display	Digital display to 1°C, adjustable in 1°C increment
	Alarm set points	Preset at -30°C for under temperature and -5°C for over temperature
Refrigerator	Control Panel	Touch-pad icon-based controls and display for easy operation Key Operated power ON/OFF switch.
	Alarms	Audio and visual alarms with battery back-up system Supplied with remote alarm contact as standard. Automatic reset and ring back functions
	Alarm conditions	High-Alarm, Low-Alarm, A.C. Failure, Low Battery, Door Open and Clean Filter
	Alarm set points	Preset at 2°C for under temperature and 8°C for over temperature
Cooling method	Aircool by circulation fan	
Safety devices	Over temperature protection, built-in overload protection for compressor, ELCB	
Refrigerators	Hermetic seal type	
Refrigerant	404A, 134a (CFC Free)	
Defrost	Auto defrost by built-in heater	
Capacity	350 litres	
Shelves	4 pcs	
Chamber dimensions (mm)	610W x 480D x 715H (Refrigerator)	610W x 480D x 525H (Freezer)
Overall dimensions (mm)	720W x 740D x 2030H	
Power source	240VAC , 50Hz, 8.2A	
Weight	130 kg	

All specifications are subject to change without prior notice.