



MODEL CSC-9500



MODEL CSC-91000



MODEL CSC-91000C

**Features**

- Total elimination of noxious vapours by ventilation fan.
- Low running cost as no expensive consumables are required e.g. sampling kits, test kits or chemical filters.
- Safe to use as there is no risk of non-functioning saturated filters.
- No hassle as there is no need to have routine or scheduled samplings or checks on filters.
- Constructed of heavy gauge mild steel plate coated with anti-corrosion powder paint.
- Doors constructed of PMMA which does not support burning in case of fire.
- Excellent visual comfort with large doors which allows full view of the stored chemicals.
- Key lock security prevents access to unauthorised persons.
- Magnetic labelling tags.
- “Fan on” indicator light.
- Adjustable fan speed.
- Silent operation.
- Adjustable shelves – Four supplied for each storage compartment.
- Supplied with spillage tray(s) and 3 meters of flexible duct piping.

**Specifications**

MODEL	CSC-9500	CSC-9500C	CSC-91000	CSC-91000C
Type	Single compartment Single door	Double compartment Double-door	Double compartment Double-door	Four compartment Four-door
Capacity	100 one litre bottles and flasks	100 one litre bottles and flasks	200 one litre bottles and flasks	200 one litre bottles and flasks
Applications	For the storage of acids, bases or solvents	For the storage of acids, bases or acid-solvents or solvents-acids in 2 isolated compartments	For the storage of acids, bases or solvents in 2 isolated compartments	For the storage of acids, bases or solvents in 4 isolated compartments
Ventilation rate	Up to 275 m <sup>3</sup> /hr	Up to 275 m <sup>3</sup> /hr	Up to 550 m <sup>3</sup> /hr	Up to 550 m <sup>3</sup> /hr
Fan speed	Adjustable	Adjustable	Adjustable	Adjustable
Outlet duct diameter	4"	4"	4"	4"
Shelf dimensions (mm)	697W x 420D x 20H (4 units)	320W x 420D x 20H (8 units)	718W x 415D x 20H (8 units)	330W x 415D x 20H (16 units)
Internal dimensions of each compartment (mm)	720W x 470D x 1690H	345W x 470D x 1690H (Two compartments)	745W x 470D x 1710H (Two compartments)	355W x 470D x 1710H (Four compartments)
External dimensions (mm)	800W x 605D x 2125H	800W x 605D x 2125H	1600W x 605D x 2150H	1600W x 605D x 2150H
Spillage tray capacity	20 litres	20 litres	20 litres (2 units)	20 litres (2 units)
Electrical supply	240 VAC, 50 Hz, 0.4 A	240 VAC, 50 Hz, 0.4 A	240 VAC, 50 Hz, 0.8 A	240 VAC, 50 Hz, 0.8 A
Weight	145 kg	155 kg	325 kg	345 kg

All specifications are subject to change without prior notice.