



Model WSC-200G

**Wall Hung Cabinets:**

The wall hung cabinets are available in modular units of size 1000mm width, 340mm depth and 700mm height. The units are constructed of 1.2mm powder coated mild steel plate. There is a choice of swing door and sliding door models. All wall-hung cabinets are supplied with brackets for purpose of hanging the units onto the walls.

**Specifications**

Type	Model	Description	Width (mm)	Depth (mm)	Height (mm)
Wall Hung Cabinets	WSC-200M	Solid door	1000	340	700
	WSC-200G	Glass door	1000	340	700

All specifications are subject to change without prior notice.



Model BSC-350-2D

**Mobile Base Cabinets:**

Constructed of 1.2mm powder coated mild steel plate, these mobile cabinets are available in modular units of size 1000mm width, 450mm depth and 770mm height (without caster wheels). You have a choice of single, double, three or four drawers. Each shelf compartment is supplied with one unit of shelf.

**Specifications**

Type	Model	Description	Width (mm)	Depth (mm)	Height (mm)
Mobile Base Cabinets	BSC-350-2D-2W	Two drawers / two doors	1000	450	770
	BSC-350-2D	Two doors	1000	450	770
	BSC-350-2D-1W	One drawer / two doors	1000	450	770

All specifications are subject to change without prior notice.

**Other related products:**

- Chemical storage cabinets.
- Acid and corrosive liquid storage cabinets.
- Flammable liquid storage cabinets.
- Stainless steel tables, racks and trolleys.
- Chemical / solution storage racks.

All specifications are subject to change without prior notice.

**Option**



Fibre Glass Tray