



Stainless steel table with SS undershelf



Stainless steel table with open base

Features

- Stainless steel construction resists rust, ensuring a clean work surface for the life of these tables.
- Ideal for applications where cleanliness is critical such as food analysis, specimen preparation and clean rooms.
- Choose from two styles. Flat-topped models offer quick, easy cleaning and feature stainless steel undershelf for storage. Raised countertop-edged tables protect floors by confining spills to tabletop.
- Heliarc welded construction creates a uniform, rock-solid unit. Tops are made of 14GA 304 stainless steel with 1 7/8" turn down on all four sides. Surface is polished to a #4 mill finish. Heavy gauge steel box channels provide rigidity and strength. Legs are 1 3/4" OD, 16 gauge 304 series solid stainless steel tubing complete with 1-in, adjustable, impact resistant, white metal feet.
- Choose from models with 18 gauge 304 series stainless steel undershelf for extra storage or open base with more unobstructed leg room.

Stainless steel table with SS undershelf

Catalogue No	Dimensions (in)
FT-SST-233D	36W x 24D x 36H
FT-SST-243D	48W x 24D x 36H
FT-SST-263D	72W x 24D x 36H
FT-SST-283D	96W x 24D x 36H
FT-SST-333D	36W x 36D x 36H
FT-SST-343D	48W x 36D x 36H
FT-SST-363D	72W x 36D x 36H
FT-SST-383D	90W x 36D x 36H

Stainless steel table with open base

Catalogue No	Dimensions (in)
SST-233D	36W x 24D x 36H
SST-243D	48W x 24D x 36H
SST-263D	72W x 24D x 36H
SST-283D	96W x 24D x 36H
SST-333D	36W x 36D x 36H
SST-343D	48W x 36D x 36H
SST-363D	72W x 36D x 36H
SST-383D	90W x 36D x 36H

All specifications are subject to change without prior notice.

Option: Raised countertop-edged table