



Model DW-FW110 / DW-FW251 / DW-FW351

Features

- Chest type
- Temperature range from -10°C to -40°C, can be set freely
- Imported microprocessor-based temperature controller
- LCD display
- Exterior body is painted steel board, interior is stainless steel (DW-FW110 is aluminium inside)
- Door with safety lock
- Balance hinge ease door opening
- Audible/visual alarm for high or low temperature
- 90 mm extra thick insulation of polyurethane
- Imported Danfoss (1 unit) compressor and EBM fan motor
- 4 casters provided for easy mobility
- Refrigerant : R134a, CFC free
- Power supply : 220-240VAC, 50 Hz, 1 phase
- Certificates : ISO9001, ISO14001, ISO13485
- Optional : basket

Specifications

Model	DW-FW110	DW-FW251	DW-FW351
Capacity	110 litres	251 litres	351 litres
Temperature	-10°C to -40°C		
Input power	260 watt	265 watt	300 watt
External size (mm)	583W x 619D x 870H	730W x 1204D x 968H	730W x 1554D x 968H
Internal size (mm)	410W x 410D x 654H	430W x 1050D x 556H	430W x 1300D x 627H
Weight (Net/gross)	34/38 kg	86/90 kg	102/110 kg

All specifications are subject to change without prior notice.