



Model DW-FW150 / DW-FW270 / DW-FW360

Features

- Chest type
- Temperature range from -20°C to -40°C, can be set freely
- Microprocessor-based temperature controller
- LED temperature display
- Exterior body is painted steel board, interior is stainless steel
- Door with safety lock
- Audible/visual alarm for high or low temperature, system failure alarm
- 110 mm thick foaming insulation
- 4 casters provided for easy mobility
- Power supply : 220-240VAC, 50 Hz, 1 phase
- Certificates : ISO9001, ISO13485

Specifications

Model	DW-FW150	DW-FW270	DW-FW360
Capacity	150 litres	270 litres	360 litres
Temperature	-20°C to -40°C		
Input power	288 watt	315 watt	360 watt
External size (mm)	811W x 775D x 929H	1245W x 775D x 929H	1534W x 775D x 929H
Internal size (mm)	585W x 465D x 651H	1019W x 465D x 651H	1308W x 465D x 651H
Weight (Net/gross)	60/67 kg	77/85 kg	98/106 kg

All specifications are subject to change without prior notice.