



Model DW-FL270

Features

- Vertical cabinet (upright)
- Operating temperature ranging from -20°C to -40°C can be set freely
- Microprocessor based temperature controller
- Digital temperature display
- Door with safety lock
- Audible / visual alarm : High or low temperature alarm, system failure alarm
- 4 casters provided for easy mobility
- Certificates : ISO9001, ISO13485

Specifications

Model	DW-FL90 (CE)	DW-FL270 (CE)	DW-FL450 (CE)
CE MARK	Yes	Yes	Yes
Capacity	90 litres	270 litres	450 litres
Drawers/shelves (fixed position)	5 drawers	7 drawers	12 drawers
Exterior construction	Powder coated steel		
Cabinet construction	ABS interior		
Insulation	90 mm foaming insulation		
External size (mm)	544W x 623D x 842H	600W x 620D x 1750H	655W x 780D x 1845H
Internal size (mm)	390W x 410D x 539H	410W x 414D x 1112H	530W x 626D x 1356H
Input power	300 watt	350 watt	650 watt
Electrical	220 – 240 VAC, 50 Hz, single phase		
Weight	45/50 kg	65/70 kg	120/125 kg

All specifications are subject to change without prior notice.