



Model CS-304

Chemstore Flammable Liquid Storage Cabinets

The Chemstore Flammable Liquid Storage Cabinets are designed for the safety storage of flammable liquid containers. These highly visible and durable yellow powder-coated cabinets feature a “HIGHLY FLAMMABLE – KEEP FIRE AWAY” warning on the front doors. These storage cabinets are constructed of double-walled, 18-gauge steel with 1½-inch (3.8 cm) insulating air barrier. Features include piano hinges, three-point lock with handle, dual vents with flash arrestors, 2-inch (5 cm) raised, leakproof door sill to contain spills and four levelling legs. Four brackets support each adjustable shelf.



Specifications

MODEL	CS-304
Capacity	4 gallons
No. of shelves supplied	1 unit
Shelf dimensions (mm)	347W x 347D x 20H
Interior dimensions (mm)	370W x 368D x 546H
Exterior dimensions (mm)	436W x 436D x 630H
Weight	35 kg

All specifications are subject to change without prior notice.

MANUFACTURED IN ACCORDANCE WITH OSHA REGULATIONS AND NFPA CODE 30 SPECIFICATIONS

VENT BOX
MODEL VBX-4



Model VBX-4

Features

- Stainless steel construction
- Each vent box can connect to 4 safety cabinets.
- 3 fan speed selections
- Fan indicator light, fan failure indicator alarm
- Choice of acid, organic, formal dehyde,filter, tyre.

Specifications

MODEL	VBX-4
Dimension of vent box (mm)	12 gallons
Dimension of exhaust part (mm)	1 unit
Dimension of inlet parts (mm) (4 pods)	520W x 380D x 35H
Fan	530W x 380D x 855H
Electrical	240 VAC, 50 Hz, 2.0A

All specifications are subject to change without prior notice.