



Model ASCPP-60

Description of item

- Freestanding single or double door models made of HDPE for storing acids, bases and corrosives.
- Cabinet completely made of HDPE, no metal parts.
- One or two manual closed doors, each door can open nearly 180° degree for easy access.
- Can be padlocked to prevent unauthorized access.
- Counter height top doubles as a convenient work surface.
- Optional cabinet ventilation opening to prevent vapor accumulation in the cabinet.

Specifications

MODEL	ASCPP-60	ASCPP-178
Capacity (Quantity x Volume)	8 bottles x 2.5 L	30 bottles x 2.5 L
No. of shelves supplied	1 unit	1 unit
Interior dimensions (mm)	352 (W) x 352 (D) x 481 (H)	790 (W) x 470 (D) x 480 (H)
Exterior dimensions (mm)	372 (W) x 382 (D) x 541 (H)	810 (W) x 500 (D) x 540 (H)

All specifications are subject to change without prior notice.